

# BIOENDO

## Lyophilized Amebocyte Lysate Products and Service



**BIOENDO**

Endotoxin  
and  
Beta-glucan  
Detection  
Expert

Dedicated to  
Endotoxin & Beta-glucan Assay  
Since the Year of 1978

**Xiamen Bioendo Technology Co., Ltd**

Xiamen Bioendo Technology Co., Ltd., founded in 1978, is the expert in the field of endotoxins and Fungi (1,3)- $\beta$ -D-glucan testing. Bioendo works on promoting application and development of Lyophilized Amebocyte Lysate, and dedicates to researching, developing, producing, and marketing Amebocyte Lysate for four decades. Our products are registered at CFDA (China Food and Drug Administration) since 1988. We participate in preparing national standard Amebocyte Lysate and Reference Standard Endotoxin and standardizing Control Standard Endotoxin for the Government.



### Gel Clot Lyophilized Amebocyte Lysate Single Test in Ampoule GS44 Series

Gel Clot Lyophilized Amebocyte Lysate Single Test in Glass ampoule endotoxin-specific Amebocyte Lysate which includes beta-glucan inhibitor in the formulation and will not react to beta-glucan. Furthermore, sample could be added to the glass ampoule directly. This means you don't need to reconstitute Amebocyte Lysate at first, and could decide how many tests use each time to avoid waste.

### Rapid Gel Clot single test kit RG Series

Rapid Gel Clot Endotoxin Assay Kit, the featured application in the field of testing water or dialysate. Easy and convenient operation, rapidly detect the endotoxins in the field of dialysis, especially it is used for the related medical device to quantify the endotoxins level. The kit is a kind of rapid solution for the medical device to test the endotoxins, no need to configure the sophisticated instrument in the testing procedure.




### Rapid Gel Clot 10 Samples Kit RG Series

Rapid Gel Clot Endotoxin Assay Kit is designed to rapidly specify endotoxin level in water and dialysate. For 10 samples kit, results could be read within 30 minutes. The 10 sample test kit is a kind of rapid solution for the medical device to test the endotoxins which provide a reasonable packing size of multi samples, specialize for the endotoxin testing in water and dialysate or related medical devices. ( Note the kit is used for the endotoxin testing in the products of medical device.)

### Gel Clot Lyophilized Amebocyte Lysate Multi-test Vial G52

Gel Clot method Lyophilized Amebocyte Lysate Multi-test Vial G52 is especially suitable for endotoxin detection in large amounts. Bioendo's lysate reagent features strong resistance to interference and strong improvement for heat stability.



### Gel Clot Lyophilized Amebocyte Lysate Multi-test Vial G17

Gel Clot Method Lyophilized Amebocyte Lysate Multi-test Vial G17 is especially suitable for endotoxin detection in large amounts. Bioendo Lysate reagent features are strong resistance to interference and strong improvement for heat stability.




### LAL Reagent Water (Water for Bacterial Endotoxins Test)

We offer various of low endotoxin, endotoxin free accessories, includes BET water (LAL Reagent Water / TAL Reagent Water), pyrogen-free test tubes, pyrogen-free pipette tips, pyrogen-free microplates, pyrogen free resevirs and so on. High quality pyrogen free products insure the success of your experiments.



### Control Standard Endotoxin (CSE)



### Pyrogen-free Pipette tips and Consumables



### Endotoxin-free Glass Test Tubes



### Pyrogen-free Centrifuge Tube

*Please contact us for more*